

Dragon Medical 360 Direct – External Access

Introduction

Dragon Medical 360 Direct is now available for external use for current Dragon Direct Users. Please Contact your local IS department for training and access requests.

System Requirements

- Windows 7 compatible PC (x64 or x86)
- High-speed internet connection (minimum 10 mbps)
- Internet Explorer 9 or above
- Citrix Receiver version 4.2 (also known as 14.2.100)
- SpeechMagic Citrix Extension Server x86 - 17.5.10.1
- PowerMic Citrix Server & Virtual Desktop Extension (32-bit) - 4.2.19.1

Hardware Requirements


Cerner and Dragon Medical 360 Direct have only been certified by UHS to work with the Nuance PowerMic II Speech Recognition Microphone which can be purchased at www.amazon.com



Installation instructions:

1. Login to Cerner and complete installs for Steps 2 & 3 using the links in the screenshot below
2. Download and install Citrix Receiver v4.2 then reboot your computer
3. Download and install the PowerMic and Speech Magic drivers then reboot your computer

Important: If your computer has a previous version of Citrix installed, you must uninstall it before installing a newer version. Please download the Citrix Cleanup Utility and run it to uninstall all previous versions of Citrix before proceeding.



[Download Receiver for Windows](#) | [Download Receiver for Mac](#) | [Download Citrix Cleanup Utility](#) | [Download Dragon Drivers](#)

*** For additional Support, please contact your normal IS support group ***